

1. Communication protocol(V5.0A)

1.1 Communication settings(for Comm port):

Data format:

1 Start bit + 8 data bits + 1 stop bit, odd;

Baudrate :

4800 baud

1.2 Real time data sent to PC from Pulse Oximeter:

DATA: 5 bytes in 1 package, 60 packages/second, bit 7 stand for synchronization。

byte	bit	content
1	0~3	Signal strength for pulsate(0~8)
	4	1= searching too long, 0=OK
	5	1=dropping of SpO2, 0=OK
	6	1=beep flag
	7	Synchronization,always be 1
2	0~6	pulse waveform data
	7	synchronization, always be 0
3	0~3	bar graph (stand for pulsate case)
	4	1=probe error, 0=OK
	5	1= searching, 0=OK
	6	bit 7 for Pulse Rate
	7	synchronization, always be 0
4	0~6	bit 0~bit 6 for Pulse Rate
	7	synchronization, always be 0
5	0~6	bit 0~bit 6 for SpO2
	7	synchronization, always be 0